# StarPhoresis

## **2-Gel Mini Vertical Electrophoresis Systems**

N2510-1010, N2516-1410, N2520-2010

## **2-Gel Wide Mini Vertical Electrophoresis System**

N2520-1010

User Manual • Manuel d'utilisation • Bedienungsanleitung • Manuale Utilizzatore

<b>English</b>	2 – 14
<b>Français</b>	15 – 27
<b>Deutsch</b>	28 – 41
<b>Italiano</b>	42 – 55

## ENGLISH

### StarPhoresis 2-Gel Mini Vertical Systems

Catalogue Nos. N2510-1010, N2516-1410, N2520-2010

### StarPhoresis 2-Gel Wide Mini Vertical System

Catalogue Nos. N2520-1010

**Please read the entire User Manual before operating this unit. These units are capable of delivering potentially lethal voltage when connected to a power supply and should be operated only by qualified and technically trained personnel.**


### Warranty

Please check that the unit has been received complete and undamaged. Refer to the illustration below and to *TABLE A* to check that all components are present. Save all packaging and documents until you have thoroughly inspected your shipment; if you find that your order is incorrect or damaged please call STARLAB for instructions for its return.

This warranty is valid for 36 months, but only if the product has been used and cared for according to this User Manual. No liability is accepted for loss or damage arising from incorrect use. The liability is limited to the repair or replacement of the unit, or a refund of the purchase price, at STARLAB's option. STARLAB is not liable for any consequential damages. We refer to our general terms for sale.

STARLAB reserves the right to alter the specifications of the Electrophoresis Systems without prior notice. This will enable us to implement improvements as soon as they are available.

### StarPhoresis 2-Gel Mini and Wide Mini Vertical Electrophoresis System components

**TABLE A:** System components included

	N2510-1010	N2516-1410	N2520-2010	N2520-1010
<b>Safety Lid</b> (with integral power cords)	1	1	1	1
<b>Upper Buffer Chamber</b> (with colour-coded electrodes, gold plated banana plugs and water cooling system)	1	1	1	1
<b>Lower Buffer Chamber</b>	1	1	1	1
<b>Casting Base</b>	1	1	1	1
<b>Blank Glass Plates (x4)</b> (W x L x Thickness)	10 x 10 x 0.24cm	16 x 14 x 0.32cm	20 x 20 x 0.32cm	20 x 10 x 0.32cm
<b>Notched Glass Plates</b> (W x L x Thickness cm)	10 x 10 x 0.24cm	16 x 14 x 0.32cm	20 x 20 x 0.32cm	20 x 10 x 0.32cm
<b>Blocking Plate</b> for running one gel	1	1	1	1
<b>Side Spacers</b>	0.8mm thick (pack of 4)	1.5mm thick (pack of 8)	1.5mm thick (pack of 8)	0.8mm thick (pack of 8)
<b>Teflon® Combs (2 of each)</b>	10 wells, 0.8mm thick	15 & 20 wells, 1.5mm thick	15 & 20 wells, 1.5mm thick	15 & 20 wells, 0.8mm thick
<b>Replacement gaskets (x2)</b> for upper buffer chamber	1	1	1	1
<b>Replacement gaskets (x2)</b> for casting base	1	1	1	1

**TABLE OF CONTENTS**

I.	ENVIRONMENTAL CONDITIONS FOR USE .....	4
II.	SAFETY PRECAUTIONS .....	4
III.	INTRODUCTION .....	4
IV.	GENERAL CARE AND CLEANING .....	5
V.	SPECIFICATIONS .....	6
VI.	UNIT ASSEMBLY FOR CASTING AND RUNNING GELS .....	6
VII.	OPERATION .....	7
VIII.	TROUBLESHOOTING .....	10
IX.	COMB SPECIFICATIONS .....	12
X.	ACCESSORIES .....	14

## I. ENVIRONMENTAL CONDITIONS FOR USE

- This unit is intended for indoor use only.
- This unit can be operated safely at an altitude of 2,000m.
- The normal operating temperature range is between 4°C and 65°C.
- Maximum relative humidity 80% for temperatures up to 31°C, decreasing linearly to 50%, relative humidity at 40°C.

## II. SAFETY PRECAUTIONS

Please read the User Manual carefully before using the STARLAB Vertical Electrophoresis Unit. This manual contains important operating and safety information. STARLAB's electrophoresis units are designed to perform flawlessly for years in the most demanding laboratories. Please take the time to read this manual so you understand the safety and operating instructions to ensure the successful use of the unit. Alterations could cause serious injury to the user or the system.

Power to the unit is supplied by an external power supply. The power supply must meet safety standards for IEC 1010-1 regulations and must be ground isolated and incorporate a no-load detecting circuit.

Power is supplied to the gel through the cover of the system providing a safety interlock to the user. Users should not attempt to operate this unit without the safety inter-locked cover in place.

- To avoid the risk of personal shock always disconnect the unit from the power supply before:  
a) removing cover , b) moving the unit or, c) adding running buffer.
- Running conditions should not exceed the maximum operating voltage or current.
- DO NOT fill the buffer chamber with running buffer above the maximum fill line.
- Use this apparatus only for its intended purpose and as described in this manual.
- Do not use the product if the power cords are damaged or if its surfaces are cracked.
- Do not operate electrophoresis units in metal trays.

**CAUTION:** During electrophoresis very low quantities of various gases are produced at the electrodes. The type of gas produced depends on the composition of the buffer employed. To disperse these gases make sure that the unit is run in a well-ventilated area.

## III. INTRODUCTION

Thank you for purchasing a StarPhoresis 2-Gel Mini or Wide Mini Vertical Electrophoresis System. Our vertical systems allow for fine resolution of protein or nucleic acid fragments on one or two acrylamide gels (PAGE). PAGE separation offers the superior resolution necessary to separate native or denatured proteins and nucleic acids in applications such as Single-Stranded-Confirmation-Polymorphism (SSCP) or dinucleotide repeat analysis using western blotting, and also for automated protein sequencing analysis. For individual sample volumes we offer a broad range of different combs with thicknesses of 0.8mm and 1.5mm.

## Outstanding features ensure trouble-free use

- Robust, acrylic construction stands up to daily usage without breakage, warping or leakage.
- All of STARLAB's Vertical Electrophoresis Systems are hand-made to provide unmatched durability for years of daily use.
- Systems meet or exceed all safety requirements for the IEC 1010-1 Standards.
- Rugged, spring-loaded clamp mechanism, alignment pins and hollow gaskets guarantee reliable leak-proof gel installation.
- Precision glass plates provide exceptional flatness and finished edges to ensure uniform separation.
- The casting base enables casting directly on the upper buffer chamber which removes the need to move gels once polymerised.

## Intelligent design results in exceptional resolution

- Electrode configuration assures uniform field, straight lanes and rapid runs, thus saving time and improving data generation rate.
- Integrated water cooling system prevents band distortion.

## IV. GENERAL CARE AND CLEANING

**WARNING!** Acrylic is not resistant to aromatic or halogenated hydrocarbons, ketones or esters. Organic solvents cause acrylic to "craze" or crack. Do not use ethanol or other organic solvents to clean the unit.

- Do not autoclave, bake or microwave the unit.
- Before using, clean and dry the unit. Clean with DISTILLED WATER ONLY; dry parts with clean, lint-free lab wipes or, preferably, air dry. Use care when cleaning or drying the unit near the platinum wire. The connectors should be clean and dry before usage or storage.
- Do not use abrasive creams or scourers.
- Do not use cleaning brushes in the electrode area.
- A thorough rinse with distilled water is all that is generally required to clean the unit after use. A mild detergent may also be used. Acrylic can also be exposed to a mild bleach solution (10:1). In addition RNase removal products are also safe for acrylic.

## Gel Plate Preparation

Clean the glass plates, spacers and combs in a mild laboratory detergent. DO NOT use abrasive creams or scourers. If a particularly clean finish is required (eg for silver-stained gels) glass plates can be soaked in chromic acid overnight, rinsed with water, then wiped successively with ethanol, acetone then ethanol again. DO NOT ALLOW organic solvents or chromic acid to come into contact with the acrylic components of your vertical system.

## Combs and Spacers

Combs and spacers for vertical systems are available in thicknesses of 0.8mm and 1.5mm. Comb spines are acrylic and the teeth are polycarbonate. Please see Chapter 9, *Comb Specifications*, for details of the combs available for each device. Combs are made of high quality Teflon®. Combs designated 'Microtitre' are compatible for use with multi-channel pipettes.


## V. SPECIFICATIONS

**TABLE B:** General specifications


	N2510-1010	N2516-1410	N2520-2010	N2520-1010
Gel dimensions (w x l)	10 x 10cm	16 x 14cm	20 x 20cm	20 x 10cm
Upper buffer chamber capacity	170ml	400ml	600ml	400ml
Lower buffer chamber capacity	~ 240ml	~ 300ml	~ 800ml	~ 450ml
Total running buffer	~ 450ml	~ 650ml	~ 1250ml	~ 750ml
Total buffer capacity	~ 450ml	~ 650ml	~ 1250ml	1100–1300ml
Current (mAmps), constant	15–35mA/Gel	15–50mA/Gel	15–75mA/Gel	30–45mA/Gel
Max. voltage (volts)	600V	600V	600V	600V
Time requirements	30–90 minutes	60–120 minutes	60–180 minutes	30–90 minutes
Sample capacity	24	48	50	72

## VI. UNIT ASSEMBLY FOR CASTING AND RUNNING GELS

1. Turn over the casting base so that the four acrylic positioning squares are on the top.
2. Place the upper buffer chamber on the casting base. The precision-machined upper buffer chamber should fit snugly over the four squares.


3. Loosen the wing nuts and slide the clamp bars outward. Please note the nylon alignment pins. These pins assure that the glass plates are properly placed over the upper buffer chamber gasket, while the acrylic squares on the casting base assure that the glass is positioned evenly and at the proper place for optimum sealing on the casting base gaskets.
4. Gel Plate assembly: On a clean, level bench position the two side spacers flush with the edges of the blank glass plate and overlay the notched plate.


Place the gel plate assemblies (cassettes) with the notched glass plate innermost against the gasket of the upper buffer chamber. The assemblies should rest on the positioning squares of the casting base and between the nylon alignment pins.

5. To clamp the assembly to the upper buffer chamber slide the clamp bars towards the middle. Tighten the wing nuts until a seal is formed between the gasket and the glass. The hollow gasket allows for a superb seal without over-tightening the wing nuts. Over-tightening may cause the glass to break.
6. If casting and running two gels, repeat steps 4 to 6 for the second gel plate assembly. If you wish to run only one gel, secure the blocking plate to the second side. Note: a combination of two gel assemblies or one gel assembly and the blocking plate are necessary to form the walls of the upper buffer chamber.

7. Lift the assembled upper buffer chamber and turn the casting base over. Turn the cams so that the handle is pointing up and pull out. Place the upper buffer chamber on the gaskets. Note: A protective plastic film is left on the gaskets for shipping. A piece of clear tape has been fastened to the end of the film to assist the user in removing it
8. Insert the cam pins and simultaneously turn the handles one half turn to tighten the assembly down onto the gasket base. Once the Upper Buffer Chamber assembly has been secured onto the casting base an initial leak test using a small amount of water is recommended – add 2ml–3ml of water and let stand for 2 minutes. If no leakage is visible, empty water and proceed with gel casting.

## VII. OPERATION

### Casting the Gel

1. To ensure reproducibility and uniform polyacrylamide cross-linking we recommend de-ionizing, de-gassing and filtration of acrylamide gel solutions prior to use. Acrylamide solutions should be stored in a cool, dark environment such as a refrigerator and allowed to reach room temperature prior to pouring. Avoid exposure to heat and sunlight. Polymerization conditions should be adjusted to effect polymerization within about 5–10 minutes. Test a small volume in a vial prior to pouring the gel. As a rough guide, 100ml of degassed 6% acrylamide gel will set in about 5 minutes at room temperature when gently mixed with 450 $\mu$ l of freshly prepared 10% (w/v) Ammonium persulfate plus 200 $\mu$ l TEMED. The setting time increases to about 10 minutes if the TEMED volume is reduced to 100 $\mu$ l and to approximately 15 minutes with 75 $\mu$ l. The amount of catalysts may need to be reduced under warm conditions. Do not pour under direct sunlight.

Prepare the appropriate volume of acrylamide gel solution using Table B below as a guide. These volumes have been calculated using the glass and spacers provided by STARLAB, and subtracting the volume of the spacers and the notch. The volumes are approximate.

**TABLE C:** Approximate gel solution volumes

	N2510-1010	N2516-1410	N2520-2010	N2520-1010
Gel solution volume, 0.8mm spacer thickness	7.5ml	13.5ml	24.6ml	15ml
Gel solution volume, 1.5mm spacer thickness	15ml	27ml	49.1ml	30ml

2. Run the acrylamide gel solution mix slowly down the inside edge of the gel cassette. **Avoid aeration.** Place the comb in the gel plate assembly.
- If a stacking gel is to be used, carefully overlay the gel solution to a depth of 3–5mm with 1 x gel buffer or water-saturated butanol. Following polymerization of the separating gel, pour off the overlay layer (rinse off butanol with electrophoresis gel buffer) and pour a stacking gel if required. Insert the comb ensuring bubbles are not trapped around comb teeth. Once the stacking gel has polymerised use the gel immediately or store wrapped in a damp paper towel and plastic film at 4C°. Wait a minimum of 15 minutes for the gel to polymerise. Repeat process as required.
3. Release the cams and pull away from the upper buffer chamber and gels. Wash off any residual acrylamide. Place the upper buffer chamber into the lower buffer chamber. Stainless steel pins are located on the lower sides of the upper buffer chamber that slide into the precision machined clear sides of the lower buffer chamber to set it in place.
4. Add the appropriate volume of running buffer to the upper buffer chamber (Table B gives approximate volumes), making sure the running buffer is 3mm below the top of the blank glass, ensuring sufficient contact with the top of the gel surface. Be sure that the running buffer is not leaking from the upper buffer chamber to the lower buffer chamber. If buffer is leaking you will need to drain the UBC and reset the gel cassettes.
- NOTE: When running only one gel a blocking plate is required on the other side of the unit to retain the top buffer level.

**Gel and buffer volumes**

Some guidelines for general operating conditions are given in Table B, but conditions vary according to the number of gels, their composition, length, and cross sectional area. The current requirement will increase in proportion to the number of gels or gel thickness providing that the voltage is not limiting, eg. two gels require twice the current of one, but the same voltage. Longer gels require proportionally higher voltages. By increasing the gel concentration the electrical resistance is increased and the rate of migration decreases. Higher voltages can be applied but be careful not to overheat the gel. The conductivity of non-dissociating buffer systems gels vary enormously and conditions must be determined empirically.

The run conditions are to be taken as a guideline only and apply to SDS Tris-glycine gels. If the plates become hot increase the water flow rates within the recommended limits or reduce the power settings. STARLAB offers a broad range of power supplies for a number of electrophoresis separations.

**Sample loading**

- If a native gel is being used, pre-electrophorese the gel for 15–40 minutes prior to loading samples. SDS gels do not need this step.
- Centrifuge samples at 12,000xg for 5 minutes. If this step is omitted samples may streak during electrophoresis.
- Carefully remove the sample comb and immediately flush the wells with electrophoresis buffer using a syringe.

- Load the samples using a gel loading pipette tip (see STARLAB's range of gel-loading tips).  
See Appendix IX. *Comb/Spacer Specification* for approximate well volumes. If possible, avoid taking liquid from the pellet area at the bottom of the tube. During sample loading the pipette tip should be 1–2mm above the bottom of the well to minimise dilution of the sample and to keep the sample as a tight layer.
- Fill unused wells with the equivalent volume of sample buffer to maintain uniform electrical resistance across the gel.
- Add buffer to the lower buffer chamber to approximately 2–3mm above the base of the gel using the Fill Line as a guide. The bottom end of the gel assembly should be in contact with the running buffer (see also VII, Gel Casting, point 4).
- Set the safety lid onto the unit so that the power cords are connected in the proper position (red to red, black to black)

## **Integrated water cooling system**

All StarPhoresis 2-Gel Electrophoresis Systems have an integrated water cooling system. Sometimes cooling is needed when running gels at a higher current or when the bioactivity of an enzyme has to be preserved. Heating of the gel can cause smearing and other problems with the resolution of protein bands. This is particularly pronounced on larger gels. We recommend running coolant or water through the cooling core in the upper buffer chamber. When ramping up voltage or current, consider at least tap water cooling.

### **To utilise Cooling System**

- Attach a separate piece of 3/8" ID clear flexible lab tubing to each hose barb on the upper buffer chamber marked as 'in' and 'out'
- Attach the tubing from the cathode (black) side of the unit, marked 'in' to either a cold-water tap or a re-circulator/chiller.
- Water flow should not exceed 2 litres per minute at 30 psi.
- Attach the drain tube on the anode (red) side marked as 'ou' to the re-circulator chiller or put into the sink drain.
- Turn on the water. Once the water has started to circulate through the system, connect the power cords to the power supply

## **Starting and ending the run**

### **Starting**

- Connect the chamber to the power supply and connect the power supply to the main electrical source. Turn all settings to zero before turning on the main source of electricity. Adjust the controls to the desired settings. Follow manufacturer's instructions.

### **Ending**

- Turn power supply settings to zero, turn off the main electrical source and disconnect the power cords. Turn off the water (if using cooling system).
- Remove the lid by pushing on the acrylic alignment pins protruding through the top of the lid with your thumbs. Slide and lift the upper buffer chamber out of the lower buffer chamber and drain buffer chambers separately.

- Loosen wing nuts and slide clamp bars outward to remove gel cassettes. It is not necessary to remove the clamp bars from the upper buffer chamber to remove the gel cassette.
- After the gel cassette(s) has been removed the gel(s) are ready for staining and blotting. Separate the plates with a strong broad blade. When using notched glass plates DO NOT pry them apart at the notches. Spread the load over a wide area.
- Rinse the chambers with distilled water then dry the electrode connectors with tissue. Ensure that the connectors are clean and dry before usage or storage.

## VIII. TROUBLESHOOTING

Many factors may affect the quality of gel preparations. For example, preparation of gel and sample buffers, gel casting and tank assembly and/or run conditions. Reading and following the instructions in this user manual can solve most problems. Some of those most commonly experienced problems, along with suggestions for solving them, are listed below:

### Problem: **ACRYLAMIDE SOLUTION LEAKS DURING CASTING**

- Ensure that the sealing surfaces of the glass plates and spacers are clean
- Ensure that each plate is free of chips
- Ensure that the wing-knobs on the upper buffer chamber are tightened (use care not to over-tighten)
- Ensure that the glass plates and spacers have been set in place using the positioning squares on the 'Flip-Side' of the gel casting base.
- Ensure that the cams in the casting base have been turned equally to tighten down the upper buffer chamber on the gaskets

### Problem: **BUBBLES DO NOT APPEAR ON THE ELECTRODES**

- Check to see if the power supply is operating properly

### Problem: **GELS FAIL TO POLYMERISE**

- May be caused by low temperatures, oxygen, insufficient/degraded catalyst, or low acrylamide concentrations

### Problem: **RUN TAKES LONGER THAN USUAL**

- Buffers may be too concentrated or at the wrong pH. Gel concentration may be too high. Check Buffer Recipe and try again. See if voltage produced by the current you are running at is the same. If it differs significantly, your buffer may not have been made up correctly.
- Upper Buffer Chamber may be leaking buffer: Make sure the gel assembly is seated firmly against the gasket. Remove gasket, wash in warm water to remove excess salts, and place the gasket back in the groove.
- Running at too low a current: Use running conditions given in this manual. When running at constant current, the current value listed is per gel.

### Problem: **RUNNING TOO FAST**

- Check buffer recipe; remake and try again. If voltage is lower than usual when running at constant current, the buffer is probably too dilute.
- Voltage or current may be set too high: turn down current setting.

Problem: **“SMILING” OF DYE FRONT**

- The centre of the gel is running hotter than at the edges: use coolant or cold tap water in cooling core and/or turn down current setting.

Problem: **VERTICAL STREAKING**

- Excessive sample or particles in the sample: either dilute sample or reduce voltage. Centrifuge samples to remove particulate contamination.
- Sample has precipitated: Centrifuge sample before adding sample buffer or use a lower percentage acrylamide gel.

Problem: **BANDS SPREAD Laterally**

- Diffusion of sample: Make sure the samples are loaded quickly and the power is applied as soon as possible after loading.
- Diffusion of sample during the run in the stacking gel: Increase percentage of stacking gel or increase current by 25% when stacking.
- Lower ionic strength of the sample: Match the ionic strength of the sample with that of the gel.

Problem: **BANDS ARE NARROWER THAN THE SAMPLE WELLS**

- Ionic strength of the sample is higher than that of the gel: De-salt the sample or use sample buffer of the same strength as the gel.

Problem: **DISTORTED SAMPLE WELLS**

- Incomplete polymerization produces poorly defined wells: De-gas gel solution prior to casting and increase APS and TEMED concentrations. The comb can be wiped with TEMED just prior to casting to improve polymerisation.
- Salt concentration is too high in the sample: Dialyse sample or use desalting column.

Problem: **RESOLVING GEL IS UNEVEN AT THE TOP**

- Overlay gel carefully using water saturated n-butanol and make sure casting stand is level.

Problem: **POORLY RESOLVED BANDS**

- May be caused by too much sample for well width or gel thickness: Dilute sample. Lower volumes generally give better resolution.
- Excessively high voltages cause fast run times, but poor resolution. Sample may have degraded.

Problem: **FROWNING OF OUTSIDE LANES**

- Leakage of Buffer along the sides or along the spacers inside the gel assembly: Do not move spacers after polymerisation and make sure that the gasket is seated firmly against the glass.

Problem: **DOUBLE BANDS – “DOUBLETS”**

- Due to re-oxidation or insufficient reduction of the sample: If using a reducing agent, prepare fresh sample buffer every 30 days. Increase the concentration of 2-mercaptoethanol or dithiothreitol in the sample.

Problem: **FEWER BANDS THAN EXPECTED WITH HEAVY BAND AT DYE FRONT**

- Caused by more than one band migrating to the dye front: increase total monomer concentration (%T).
- Sample may have degraded due to incorrect storage and/or contamination.

## IX. COMB/SPACER SPECIFICATIONS

### N2510-1010 (10 x 10cm)

COMBS				
STARLAB Cat No.	Number of Teeth	Thickness (mm)	Width (mm)	Well Volume (μl)
N2510-1056	6	0.8	11.1	142
N2510-1057	8	0.8	7.7	99
N2510-1058	10	0.8	5.7	73
N2510-1059	12	0.8	4.3	55
N2510-1026	6	1.5	11.1	266
N2510-1028	8	1.5	7.7	185
N2510-1029	10	1.5	5.7	136
N2510-1032	12	1.5	4.3	103
SPACERS				
N2510-1092	0.7mm wide x 0.8mm thick, pack of 2			
N2510-1091	0.7mm wide x 1.5mm thick, pack of 2			

### N2516-1410 (16 x 14cm)

COMBS				
STARLAB Cat No.	Number of Teeth	Thickness (mm)	Width (mm)	Well Volume (μl)
N2516-1440	10	0.8	10.4	183
N2516-1467	15	0.8	6.1	107
N2516-1437	20	0.8	3.9	69
N2516-1439	24	0.8	2.9	51
N2516-1429	10	1.5	10.4	343
N2516-1465	15	1.5	6.1	201
N2516-1430	20	1.5	3.9	129
N2516-1431	24	1.5	2.9	96
N2516-1450 Preparative comb	2	1.5	119.7/4.7	3630/152
SPACERS				
N2516-1492	1.0mm wide x 0.8mm thick, pack of 2			
N2516-1491	1.0mm wide x 1.5mm thick, pack of 2			

**N2520-2010** (20 x 20cm)

<b>COMBS</b>				
<b>STARLAB Comb Cat No.</b>	<b>Number of Teeth</b>	<b>Thickness (mm)</b>	<b>Width (mm)</b>	<b>Well Volume (µl)</b>
<b>N2520-2040</b>	10	0.8	13.6	239
<b>N2520-2067</b>	15	0.8	8.2	144
<b>N2520-2037</b>	20	0.8	5.5	97
<b>N2520-2039</b>	25	0.8	3.9	69
<b>N2520-2029</b>	10	1.5	13.6	449
<b>N2520-2065</b>	15	1.5	8.2	271
<b>N2520-2030</b>	20	1.5	5.5	182
<b>N2520-2031</b>	25	1.5	3.9	129
<b>N2520-2050</b> Preparative comb	2	1.5	148.1/4.7	4885/155
<b>SPACERS</b>				
<b>N2520-2092</b>	1.0mm wide x 0.8mm thick, pack of 2			
<b>N2520-2091</b>	1.0mm wide x 1.5mm thick, pack of 2			

**N2520-1010** (20 x 10cm)

<b>COMBS</b>				
<b>STARLAB Comb Cat No.</b>	<b>Number of Teeth</b>	<b>Thickness (mm)</b>	<b>Width (mm)</b>	<b>Well Volume (µl)</b>
<b>N2520-1040</b>	10	0.8	13.6	239
<b>N2520-1067</b>	15	0.8	8.2	144
<b>N2520-1037</b>	20	0.8	5.5	97
<b>N2520-1039</b>	25	0.8	3.9	69
<b>N2520-1073</b> Microtitre Comb	18	0.8	6.5	78
<b>N2520-1074</b> Microtitre Comb	36	0.8	2.7	32
<b>N2520-1029</b>	10	1.5	13.6	449
<b>N2520-1065</b>	15	1.5	8.2	271
<b>N2520-1030</b>	20	1.5	5.5	182
<b>N2520-1031</b>	25	1.5	3.9	129
<b>N2520-1071</b> Microtitre Comb	18	1.5	6.5	156
<b>N2520-1072</b> Microtitre Comb	36	1.5	2.7	64
<b>SPACERS</b>				
<b>N2520-1092</b>	1.2mm wide x 0.8mm thick, pack of 2			
<b>N2520-1091</b>	1.2mm wide x 1.5mm thick, pack of 2			

## X. ACCESSORIES

### Accessories for N2510-1010 (10 x 10cm)

<b>N2510-1018</b>	StarPhoresis Casting Base System, 10 x 10cm
<b>N2510-1089</b>	Replacement Gaskets for Upper Buffer Chamber, Pack of 2
<b>N2510-1019</b>	Replacement Clamps, black, Pack of 2
<b>N2510-1020</b>	Wing Knobs, Pack of 4
<b>N2510-1021</b>	Notched Glass Plate, Clear, 10 x 10cm, 0.24cm thick
<b>N2510-1022</b>	Plain Glass Plate, Clear, 10 x 10cm, 0.24cm thick
<b>N2510-1023</b>	Notched Aluminium Plate, White, 10 x 10cm, 0.24cm thick

### Accessories for N2516-1410 (16 x 14cm)

<b>N2516-1418</b>	StarPhoresis Casting Base System, 16 x 14cm
<b>N2516-1489</b>	Replacement Gaskets for Upper Buffer Chamber, Pack of 2
<b>N2516-1419</b>	Replacement Clamps, black, Pack of 2
<b>N2510-1020</b>	Wing Knobs, Pack of 4
<b>N2516-1421</b>	Notched Glass Plate, Clear, 16 x 14cm, 0.32cm thick
<b>N2516-1422</b>	Plain Glass Plate, Clear, 16 x 14cm, 0.32cm thick

### Accessories for N2520-2010 (20 x 20cm)

<b>N2520-2018</b>	StarPhoresis Casting Base System, 20 x 20cm
<b>N2520-2089</b>	Replacement Gaskets for Upper Buffer Chamber, Pack of 2
<b>N2520-2019</b>	Replacement Clamps, black, Pack of 2
<b>N2510-1020</b>	Wing Knobs, Pack of 4
<b>N2520-2021</b>	Notched Glass Plate, Clear, 20 x 20cm, 0.32cm thick
<b>N2520-2022</b>	Plain Glass Plate, Clear, 20 x 20cm, 0.32cm thick

### Accessories for N2520-1010 (20 x 10cm)

<b>N2520-1018</b>	StarPhoresis Casting Base System, 20 x 10cm
<b>N2520-1089</b>	Replacement Gaskets for Upper Buffer Chamber, Pack of 2
<b>N2520-1019</b>	Replacement Clamps, black, Pack of 2
<b>N2510-1020</b>	Wing Knobs, Pack of 4
<b>N2520-1021</b>	Notched Glass Plate, Clear, 20 x 10cm, 0.24cm thick
<b>N2520-1022</b>	Plain Glass Plate, Clear, 20 x 10cm, 0.24cm thick

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**STARLAB GmbH**
**STARLAB (UK), Ltd**
**STARLAB S.r.l.**